

PRODUCT CODES



smartflow
CONTINENCE SYSTEMS





Contents

URINARY BAGS

- Bendi Bag 4
- Comfort Bag 5
- Discreet Bag 6
- Urinary Bag (Accessories) & p.grip 7

URINARY SHEATHS

- p.sure self-adhesive latex free 8
- Comfort plus latex free dry sheaths with hydrocolloid liner 8
- Comfort latex dry sheaths 9
- Comfort Extra Thin latex dry sheaths 9
- K+ICS latex dry sheaths with removable tips 10
(for Intermittent Self Catheterisation - ISC)

URINARY SHEATHS (ACCESSORIES)

- Skin Adhesives and double sided hydrocolloid adhesive tape/liner 11
- Skin preparation and protecting fluid – Preventox 11
- Protective Cloths 12
- Drip Urinal Tap 12
- Scissors 12

URINALS

- The Uribag, Uribag F and Urifem 12

NEPHROSTOMY DRAINAGE SYSTEM 13

AIDS for CATHETERISATION 14

Bendi Bag

Community prescription packs contain:
 10 drainage bags, 1 wide strap (LB.W**), 1 narrow strap (LB.N**), connectors to attach to a night bag (55.38**) and directions for use.



1300ml



700ml

Prescription FP10 code	Sterile	Outlet tap	Capacity (ml)	Inlet length	PIP Code (required by chemist)	NHS PASA code***
70.33-04	N	lever	700	direct	297-4954	n/a
70.33-12	N	lever	700	12cm	214-2503	n/a
70.33-20	N	lever	700	20cm	214-2511	n/a
70.33-35adj	N	lever	700	35cm	297-6017	n/a
70.33-35Padj	N	sliding	700	35cm	297-6009	n/a
70.33-35Sadj	Y	lever	700	35cm	297-5985	n/a
70.33-04S	Y	lever	700	direct	297-5993	n/a
70.33-12S	Y	lever	700	12cm	214-2545	GDW004
70.33-20S	Y	lever	700	20cm	214-2552	GDW005
70.47-04	N	lever	1300	direct	297-6058	n/a
70.47-12	N	lever	1300	12cm	214-2529	n/a
70.47-20	N	lever	1300	20cm	214-2537	n/a
70.47-35adj	N	lever	1300	35cm	297-6041	n/a
70.47-35Padj	N	sliding	1300	35cm	297-6033	n/a
70.47-35Sadj	Y	lever	1300	35cm	297-6025	n/a
70.47-04S	Y	lever	1300	direct	297-6066	n/a
70.47-12S	Y	lever	1300	12cm	214-2560	GDW007
70.47-20S	Y	lever	1300	20cm	214-2578	GDW013

**available separately on prescription

***NHS PASA packs do not contain straps (order separately). Each bag has a night bag connector on each outlet tap

Comfort Bag

Community prescription packs contain
 10 drainage bags, 1 wide strap (LB.W**), 1 narrow strap (LB.N**), connectors to attach
 to a night bag (55.38**) and directions for use.




Prescription FP10 code	Sterile	Outlet tap	Capacity (ml)	Inlet length	PIP Code (required by chemist)	NHS PASA code***
710.1104	N	Lever	600	direct	270-0128	n/a
710.1204S	Y	Lever	600	direct	270-0060	n/a
710.1315adj	N	Lever	600	15cm	270-0136	n/a
710.1415Sadj	Y	Lever	600	15cm	270-0078	n/a
710.1345adj	N	Lever	600	45cm	270-0144	n/a
710.1445Sadj	Y	Lever	600	45cm	270-0086	n/a
R710.1445SOadj	Y	Reverse Lever	600	45cm	297-6082	GDW025
710.2104	N	Sliding	600	direct	270-0151	n/a
710.2204S	Y	Sliding	600	direct	270-0094	n/a
710.2315adj	N	Sliding	600	15cm	270-0169	n/a
710.2415Sadj	Y	Sliding	600	15cm	270-0102	n/a
710.2345adj	N	Sliding	600	45cm	270-0177	n/a
710.2445Sadj	Y	Sliding	600	45cm	270-0110	n/a
R711.1445SOadj	Y	Reverse Lever	1000	45cm	297-6074	n/a

**available separately on prescription

***NHS PASA packs do not contain straps (order separately). Each bag has a night bag connector on each outlet tap

Discreet Bag

Discreet Bags community prescription packs contain 10 drainage bags, 1 wide strap (LB.W**), 1 narrow strap (LB.N**), connectors to attach to a night bag (55.38**) and directions for use.



Prescription FP10 Code	Sterile	Outlet tap	Capacity (ml)	Inlet length & position	PIP Code (required by chemist)	NHS PASA Code ***
70.04S	Y	lever	450	direct corner	238-9617	GDW014
70.06S	Y	lever	400	direct centre	238-9591	GDW015
70.04	N	lever	450	direct corner	238-9583	n/a
70.06	N	lever	400	direct centre	238-9575	n/a
Mini70.04S	Y	lever	275	direct corner	249-1280	n/a

**available separately on prescription

***NHS PASA packs do not contain straps (order separately). Each bag has a night bag connector on each outlet tap

Urinary Bags (Accessories)

Comfort Leg Bag Straps & other related accessories

Prescription FP10 code	Description	Width (mm)	Qty	PIP Code (required by chemist)	NHS PASA code
LB.N	Narrow leg bag strap	20	1	304-6588	n/a
LB.W	Wide leg bag strap	38	1	304-6620	n/a
n/a	Narrow leg bag strap	20	10	n/a	GDW030
n/a	Wide leg bag strap	38	10	n/a	GDW027
55.38-10	Leg bag to night bag connector	n/a	10	304-6547	n/a

p.grip catheter/tubing retaining strap

p.grip anchors urinary catheters or drainage bag tubing to a strap to prevent drag from the appliance.

p.grip is ideal to prevent drag from day or night drainage bags (leg bags or bed bags).

Prescription FP10 code	Description	Length/Width (mm)	Qty	PIP Code (required by chemist)
CS.01	Leg strap, 2 adhesive tabs, 1 Velcro tab	600 x 38	1	304-6562
CS.02	Abdominal strap, 2 adhesive tabs, 1 Velcro tab	900 x 38	1	304-6554
CS.TABS	10 x Extra adhesive tabs, 1 Velcro tab	n/a	10	304-6570

Urinary Sheaths

p.sure self adhesive latex free sheaths

Community packs contain 30 sheaths, 1 protective cloth and directions for use.

PASA packs contain 5 sheaths, 5 tape measures, 5 adhesive patient record labels and 5 directions for use.

Prescription FP10 Code	Sheath Diameter (mm)	Prescription/FP10 Qty	PIP Code (required by chemist)	NHS PASA code	NHS PASA Qty
97.18	18	30	290-2849	GDK014	5
97.20	20	30	290-2856	GDK018	5
97.22	22	30	290-2864	GDK020	5
97.24	24	30	290-2872	GDK022	5
97.26	26	30	290-2880	GDK025	5
97.28	28	30	290-2898	GDK027	5
97.30	30	30	290-2906	GDK028	5
97.32	32	30	290-2914	GDK037	5
97.35	35	30	290-2922	GDK038	5
97.37	37	30	290-2930	GDK039	5
97.40	40	30	290-2948	GDK040	5

Comfort plus latex free dry sheaths with hydrocolloid liner

Community packs contain 30 sheaths, 30 double sided hydrocolloid adhesive tapes, 1 protective cloth and directions for use.

PASA packs contain 5 sheaths, 5 hydrocolloid adhesive liners, 5 tape measures, 5 adhesive patient record labels and 5 directions for use.

Prescription FP10 Code	Sheath Diameter (mm)	Prescription/FP10 Qty	PIP Code (required by chemist)	NHS PASA code	NHS PASA Qty
56.18	18	30	249-1215	GCE072	5
56.20	20	30	290-0918	GCE063	5
56.22	22	30	290-0926	GCE064	5
56.24	24	30	290-0934	GCE065	5
56.26	26	30	290-0942	GCE066	5
56.28	28	30	290-0959	GCE067	5
56.30	30	30	290-0967	GCE068	5
56.32	32	30	290-0975	GCE069	5
56.35	35	30	290-0983	GCE070	5
56.37	37	30	290-0991	GCE073	5
56.40	40	30	290-1007	GCE074	5

Urinary Sheaths

Comfort latex dry sheaths

Community packs contain 30 sheaths, 1 protective cloth and directions for use. Dry sheaths require a separate adhesive tape similar to our product code 50.14 which is a hydrocolloid double sided adhesive liner or a liquid skin adhesive similar to our product codes 50.01 or 50.36

Prescription FP10 Code	Sheath Diameter (mm)	Prescription/FP10 Qty	PIP Code (required by chemist)
53.18	18	30	238-9468
53.20	20	30	238-9476
53.22	22	30	238-9484
53.24	24	30	238-9492
53.26	26	30	238-9500
53.28	28	30	238-9518
53.30	30	30	238-9526
53.32	32	30	238-9534
53.35	35	30	238-9542
53.37	37	30	238-9559
53.40	40	30	238-9567

Comfort Extra Thin latex dry sheaths


Community packs contain 30 sheaths and 1 protective cloth. Dry sheaths require a separate adhesive tape similar to our product code 50.14 which is a hydrocolloid double sided adhesive liner or a liquid skin adhesive similar to our product codes 50.01 or 50.36

Prescription FP10 Code	Diameter (mm)	Prescription/FP10 Qty	PIP Code (required by chemist)
53.26D	26	30	238-9625
53.30D	30	30	238-9633
53.35D	35	30	238-9641
53.40D	40	30	238-9658

Urinary Sheaths - The K+IC System

K+ICS latex dry sheaths with removable tips for Intermittent Self Catheterisation (ISC)

Community packs contain 30 sheaths, 5 black plastic connectors, 1 protective cloth and directions for use. Dry sheaths require a separate adhesive tape similar to our product code 50.14 which is a hydrocolloid double sided adhesive liner or a liquid skin adhesive similar to our product codes 50.01 or 50.36.



Prescription FP10 Code	Sheath Diameter (mm)	Qty	PIP Code (required by chemist)
100.01	Sheath expander	1	256-3633
100.05	Extra black plastic connectors	10	244-3497
103.18	18mm	30	256-3641
103.20	20mm	30	256-3658
103.22	22mm	30	256-3666
103.24	24mm	30	256-3674
103.26	26mm	30	256-3682
103.28	28mm	30	256-3690
103.30	30mm	30	256-3708
103.32	32mm	30	256-3716
103.35	35mm	30	256-3724
103.37	37mm	30	256-3732
103.40	40mm	30	256-3740



Urinary Sheaths (Accessories)

Skin Adhesives and double sided hydrocolloid adhesive tape/liner

Required for adhering dry sheaths.

Liquid skin adhesives available latex or synthetic in a tube with pipette nozzle applicator or bottle with brush in lid.

Hydrocolloid double sided adhesive tape uses a pressure sensitive adhesive and is recommended for sensitive skin.


Prescription FP10 Code	Product name & description	Applicator	Volume/Qty	PIP Code (required by chemist)
50.00	Latex pure - no skin care component	Long pipette nozzle	2 x 28g tube	273-0018
50.01	Latex original - standard strength	Long pipette nozzle	2 x 28g tube	238-9666
50.03	Latex lanolin free - for lanolin allergies	Long pipette nozzle	2 x 28g tube	272-8343
50.05	Latex 50% - 50% less skin care	Long pipette nozzle	2 x 28g tube	272-8376
50.20	Latex 2% resin - 2% extra resin	Long pipette nozzle	2 x 28g tube	272-8392
50.22	Latex 12%resin - 12% extra resin	Long pipette nozzle	2 x 28g tube	272-8400
50.36	Synthetic latex free in bottle	Brush in lid	45ml bottle	249-1165
50.36-2	Synthetic latex free in tube	Long pipette nozzle	2 x 28g tube	249-1181
50.14	Double sided hydrocolloid adhesive tape/liner	Tape/liner	30 strips	249-1090

Skin preparation and protecting fluid – Preventox

Clean and protect skin with clear breathable film.

Increases adhesive strength of any skin to appliance bond by preparing skin before applying adhesive.

Prescription FP10 Code	Product name & description	Volume/Qty	PIP Code (required by chemist)
50.50	Preventox skin protecting wipes	50 x wipes	249-1249
50.58	Preventox skin protecting fluid - roll-on applicator	50ml bottle	249-1231



Urinary Sheaths (Accessories)...continued



Order Code	Product name & description	Available on prescription/FP10	Unit of sale	PIP Code (required by chemist)
52.70-18	Protective Cloths 18mm	N	30	306-1876
52.70-20	Protective Cloths 20mm	N	30	306-1900
52.70-22	Protective Cloths 22mm	N	30	306-1918
52.70-24	Protective Cloths 24mm	N	30	306-1926
52.70-26	Protective Cloths 26mm	N	30	306-1934
52.70-28	Protective Cloths 28mm	N	30	306-1942
52.70-30	Protective Cloths 30mm	N	30	306-1959
52.70-32	Protective Cloths 32mm	N	30	306-1967
52.70-35	Protective Cloths 35mm	N	30	306-1975
52.70-37	Protective Cloths 37mm	N	30	306-1991
52.70-40	Protective Cloths 40mm	N	30	306-1983
57.67-3	Drip Urinal Tap	N	3	306-3534
53.50	Scissors	N	each	306-3823

Call our sales line for updates on availability on Prescription/FP10

Urinals

The Uribag and Urifem

Prescription FP10 Code	Description	Capacity	PIP Code (required by chemist)
URlbag	URlbag - Male pocket-sized urinal	1100ml	249-1264
URlbag F	URlbag - Female pocket-sized urinal	1100ml	272-1140
URlfem	URlfem - Bottle urinal for ladies lying down	1200ml	249-1256

Nephrostomy Drainage System

NephSys

For use with nephrostomy catheters only with luer lock connectors. The complete system comprises of a waist belt which supports up to two drainage bag suspenders that can be positioned anywhere on the belt. The suspenders securely attach special drainage bags with luer lock connectors. Connector to attach the outlet tap to any UK night drainage bag also supplied.



Product code	Description	Inlet	Qty	PIP Code (required by chemist)	NHS PASA code
NephSys.01S	SYSTEM: 20cm inlet with luer lock connector SMALL (600mm x 38mm) waistbelt and 2 adjustable bag suspenders	20cm	4	305-1554	n/a
NephSys.01L	SYSTEM: 30cm inlet with luer lock connector SMALL (600mm x 38mm) waistbelt and 2 adjustable bag suspenders	30cm	4	305-1570	n/a
NephSys.02S	SYSTEM: 20cm inlet with luer lock connector ADULT (900mm x 38mm) waistbelt and 2 adjustable bag suspenders	20cm	4	305-1588	GDW033
NephSys.02L	SYSTEM: 30cm inlet with luer lock connector ADULT (900mm x 38mm) waistbelt and 2 adjustable bag suspenders	30cm	4	305-1562	GDW034
NS721.1720S	BAGS: drainage bag with 20cm inlet and luer lock connector	20cm	10	305-1547	n/a
NS721.1730S	BAGS: drainage bag with 30cm inlet and luer lock connector	30cm	10	305-1539	n/a
NSBelt.01	BELT: SMALL waistbelt (600x38mm) and 2 drainage bag suspenders	n/a	1	305-1521	n/a
NSBelt.02	BELT: ADULT waistbelt (900x38mm) and 2 drainage bag suspenders	n/a	1	305-1513	n/a

Aids for Catheterisation

These products are specifically designed to aid both consumers and healthcare professionals in the process of catheterisation whether permanent urethral or intermittent self-catheterisation.



110.01 standard
110.02 large

Order Code	Product name & description	Available on prescription/FP10	Unit of sale	PIP Code (required by chemist)
110.01	Inflatable leg dividers standard (mirror separate)	N*	each	305-1653
110.02	Inflatable leg dividers large (mirror separate)	N*	each	305-1661
110.05	Metal leg dividers (mirror separate)	N	each	305-1646
110.31	Mirror and strap for inflatable leg divider	N*	each	305-1679
110.32	Mirror for metal leg divider	N	each	305-1695
110.36	Torch for metal leg divider	N	each	305-1687
110.46	Labia spreader for female catheterisation	N	each	305-1638
110.50	Cath Hand for right handed person	N	each	305-1604
110.55	Cath Hand for left handed person	N	each	305-1596
110.60	p.hold penis stretcher and holder standard	N*	each	305-1612
110.61	p.hold penis stretcher and holder large	N*	each	305-1620


**Call our sales line for updates on availability on Prescription/FP10*



110.46



110.50 (right hand)
110.55 (left hand)



110.05



110.60 standard
110.61 large





Smartflow
Manfred Sauer UK
Unit 3, io Centre, Lodge Farm, Northampton NN5 7UW
Tel: 0870 1 904 100 Fax: 0870 1 904 130
Helpline: 0870 1 904 150
email: helpline@manfredsauer.co.uk www.manfredsauer.co.uk

SMARTflow® is a registered trademark of Manfred Sauer GmbH