

Our IQCath is the perfect partner for the K+IC System

IQCATH
Made by Manfred Sauer

ISC Catheter with intelligent tip
Soft spherical tip and flexible zone to navigate bends in the urethra.

K+IC

The sheath with a removable tip!

CLINICAL DEMAND

When there is a clinical requirement to intermittently self catheterise (ISC), the KIC System brings together the security of wearing a sheath without interfering with the ISC procedure.

This has been developed for patients who are incontinent between catheterisation intervals, or who achieve continence only with medication which they cannot, or will not, tolerate because of side effects. KIC overcomes the need to change sheaths after each catheterisation.

In the KIC System, the black plastic connector can be removed and the sheath, with the expander open, can be pulled over the penis so that the penis tip is left free for disinfecting and catheterisation. The sheath can then be restored to its original position and connected to the existing drainage system. This procedure can be repeated any number of times. With the KIC System you only need one sheath per day which helps protect the skin.

| Size | Order No. Prescriptions | Description |
|------|-------------------------|---|
| | 100.01 | KIC Sheath expander x 1- can be used for several months |
| | 100.05 | KIC Connector x 10, allows you to separate and reconnect leg bag from the sheath |
| 18mm | 103.18 | KIC Sheaths*, latex dry sheath, require skin adhesive or adhesive tape to secure sheath (see Accessories part of leaflet). Reinforced aperture to allow it to return to its original shape after stretching with expander over penis. * KIC Sheath packs available on doctor's prescription contain 30 sheaths, 5 black plastic connectors (100.05), 1 protective cloth to keep pubic hair away from the adhesive area of the sheath and directions for use. |
| 20mm | 103.20 | |
| 22mm | 103.22 | |
| 24mm | 103.24 | |
| 26mm | 103.26 | |
| 28mm | 103.28 | |
| 30mm | 103.30 | |
| 32mm | 103.32 | |
| 35mm | 103.35 | |
| 37mm | 103.37 | |
| 40mm | 103.40 | |

HOT TIP 1: Use our tubing gripper strap p.grip to prevent leg or night bag dragging on sheath - ask doctor to prescribe p.grip order code: CS.01

- 18
- 20
- 22
- 24
- 26
- 28
- 30
- 32
- 35
- 37
- 40

Sheath diameter in millimetres

smartflow
CONTINENCE SYSTEMS

p.grip

SMARTflow sheaths have been developed and manufactured using a unique 'dipping' process ensuring comfort, security, reliability and easy use.

11 SIZES

Means an exact fit can be achieved as there are 11 finely graded sizes (18-40mm diameter). The better the fit the safer the adhesive contact. Call for a free measuring tape to ensure you get the exact size.

THE HIGH ELASTICITY OF MATERIAL

Allows for large amounts of urine to be passed without "blowing off" the sheath, this is due to the unique "buffer zone" part of the sheath that can actually expand to hold 50mL to 100mL of urine. This also prevents urine residue remaining in the sheath after voiding reducing the risk of irritation and infection which often happens with sheaths made of a stiffer material and "funnel type" design. In fact if the user only "dribbles" a tap can be used on the sheath outlet rather than connecting the sheath to a leg bag (Order No.CVS see Accessories section of this leaflet).

THE QUALITY OF THE MATERIAL

Enables the sheath to be made of thin, but tough, gossamer which moulds itself exactly to the penis resulting in better adhesion and ultimately comfort for the wearer. A man's "shape" changes throughout the day depending on many factors SMARTflow is designed with this in mind.

THE CONNECTING TUBE

Will connect to any UK day or night bag fitted with a standard conical universal adaptor.

QUALITY ASSURANCE

All our sheaths are individually tested before leaving the factory. Safety is therefore assured.

ACCESSORIES

SMARTflow SKIN ADHESIVES

This liquid skin adhesive gives the best adhesion between skin and sheath. Unlike the self-adhering sheath, these sheaths allow each user to determine the quantity used, the width of the adhering area and the size of the buffer zone. Our skin adhesive has been used successfully for over 30 years, offering the best solution to all types of skin.

SMARTflow Original (50.01)

Latex skin adhesive reliable and secure for most skins.



Long nozzle applicator

| Latex Skin adhesives | Order no. | Strength (Standard 1 - 6 strong) |
|-------------------------------|-----------|----------------------------------|
| No skin care components | 50.00 | 4 |
| Original - standard strength | 50.01 | 1 |
| Lanolin free | 50.03 | 2 |
| 50% less skin care components | 50.05 | 3 |
| 2% extra resin | 50.20 | 5 |
| 12% extra resin | 50.22 | 6 |

SMARTflow SYNTHETIC SKIN ADHESIVES*



Bottle with applicator brush in lid



PREVENTOX SKIN PROTECTION
Protects sensitive skin

Protective cloths (50.70 - xx pack of 30) where xx is diameter of sheath

Single wipes 50.50 (pack of 50) Fluid 50.58 (roll on applicator)

Double-sided adhesive tape (50.14 - pack of 30 tapes)



Sheath tap or Catheter valve
Tap can be used on the end of sheath instead of leg bag for light incontinence (CVS - pack of 5) see HOT TIP 3

Scissors Special tips to remove sheath collar 53.50
Free Measuring tape

Call for FREE sample

MANFRED SAUER UK Ltd, Unit 3, io Centre, Lodge Farm, Northampton NN5 7UW
Tel: 01604 588090 Fax: 01604 588091 Helpline: 01604 595696
email: helpline@manfredsauer.co.uk www.urinarysheaths.co.uk

SMARTflow is a registered trademark of Manfred Sauer UK Ltd

Other SMARTflow leaflets available include: Drainage Bags • Urinals • p.grip Catheter/Tubing Retaining Strap • Nephrostomy Drainage System

SMARTflow p.sure LATEX FREE SELF-ADHESIVE SHEATHS

p.sure offers the same features as other SMARTflow sheaths but with self-adhesion

1 Correct Sizing

TO MEASURE: put tape around the middle of the penis. If the sheath is likely to be fitted when the penis is firm it should be measured when it is firm.



7 easy to follow steps for fitting p.sure sheaths

HOT TIP 2: See the trouble shooting guide in our free tutorial booklet for any issues regarding sheath fitting

p.sure

| SIZE | ORDER No. PRESCRIPTIONS* |
|------|--------------------------|
| 18mm | 97.18 |
| 20mm | 97.20 |
| 22mm | 97.22 |
| 24mm | 97.24 |
| 26mm | 97.26 |
| 28mm | 97.28 |
| 30mm | 97.30 |
| 32mm | 97.32 |
| 35mm | 97.35 |
| 37mm | 97.37 |
| 40mm | 97.40 |

*p.sure community packs available on doctor's prescription contain 30 sheaths, 1 protective cloth to keep pubic hair away from the adhesive area of the sheath and directions for use.

HOT TIP 3: If the user "dribbles" less than 100ml, a tap can be used on the sheath outlet rather than connecting the sheath to a leg bag - Order No.CVS - see Accessories section of this leaflet.

p.sure is a registered trademark of Manfred Sauer UK Ltd.

MANFRED SAUER

Not 3, 4 or 5 but 11 sizes
SHEATHS

smartflow
CONTINENCE SYSTEMS

p.sure
K+IC

SMARTflow 'COMFORT' SHEATHS

SMARTflow COMFORT PLUS SHEATHS (Synthetic)



The most hypoallergenic latex free penile sheath solution available.

Dry latex free sheath with pressure sensitive hydrocolloid double sided adhesive tape or liner

Hydrocolloid is a material used in many wound dressings and consequently is very kind to the skin and unlikely to cause allergic reactions. The skin adhesives used in self adhesive sheaths and the separate liquid skin adhesives use a solvent to prevent the adhesive setting until exposed to air, it is this solvent which often causes skin problems that can be avoided by using this type of pressure sensitive adhesive liner.

*Comfort Plus community packs available on doctor's prescription contain 30 sheaths, 30 hydrocolloid liners, 1 protective cloth to keep pubic hair away from the adhesive area of the sheath and directions for use.

| SIZE | ORDER No. PRESCRIPTIONS* |
|------|--------------------------|
| 18mm | 56.18 |
| 20mm | 56.20 |
| 22mm | 56.22 |
| 24mm | 56.24 |
| 26mm | 56.26 |
| 28mm | 56.28 |
| 30mm | 56.30 |
| 32mm | 56.32 |
| 35mm | 56.35 |
| 37mm | 56.37 |
| 40mm | 56.40 |

SMARTflow COMFORT SHEATHS (Latex)



Dry latex sheath that requires a double sided adhesive tape or separate liquid skin adhesive to attach to the skin (see Accessories section of this leaflet)

These sheaths are made from natural rubber/latex which is one of the most elastic materials available. They are less expensive to produce and can be made "thinner". Consequently they are often more secure and comfortable to wear.

Latex allergies do occur but certainly in relation to medical devices it is often the impurity of the latex used (chemicals used during the manufacturing process which are not effectively removed) that causes these allergies rather than the latex itself. We use only the highest quality materials to reduce the chance of any allergic reactions.

| SIZE | ORDER No. | |
|------|--------------------------|--------|
| | Standard Extra thin | |
| 18mm | 53.18 | |
| 20mm | 53.20 | |
| 22mm | 53.22 | |
| 24mm | 53.24 | |
| 26mm | 53.26 | 53.26D |
| 28mm | 53.28 | |
| 30mm | 53.30 | 53.30D |
| 32mm | 53.32 | |
| 35mm | 53.35 | 53.35D |
| 37mm | 53.37 | |
| 40mm | 53.40 | 53.40D |

SMARTflow COMFORT SHEATHS (Extra Thin)

These sheaths are finer than a contraceptive sheath and literally fit like a second skin. Due to their elasticity the 4 sizes available cover a 26mm to 44mm size range. Recommended for men who experience penile retraction.